





Thermoplastics


REFERENCE GUIDE


CIVCO
Radiotherapy™

CE

<p>Type-S™</p> 	 <p>Latched 已锁定 Zaključano Zachycen západkou Lást Gesloten Lukustatud Lukittu Enclenché Arretiert Ασφαλισμένο Rögzített Chiuso</p> <p>掛かっている Fiksēts Užfiskuota Lást Zabezpieczony zatraskiem Preso Preso Inciuat So západkou Zaskočen Bloqueado Rögzítetlen Mandali kapalı</p>	 <p>Unlatched 已解锁 Otključano Ödroyeno Åben Geopend Lukustamata Auki Déverrouillé Gelöst Απασφαλισμένο Rögzítetlen Aperto</p> <p>外れている Nefiksēts Neuzfiksiskuota Ulást Niezatrasnięty Nåo preso Solto Descuiat Bez západky Odpahņjen Desbloqueado Ospårrad Mandali açık</p>
---	---	--

<p>Reloadable / Uni-frame®</p> 	<p>可重复安装 / Uni-frame® Za višekratnu upotrebu / Uni-frame® Vyměnitelný/Uni-frame® Formgivningsramme til genbrug/Uni-frame® Herlaadbaar / Uni-frame® Taaslaetav / Uni-frame® Ludelleenladattava / Uni-frame® Cadre rechargeable / Uni-frame® Wiederaufladbarer Rahmen / Uni-frame® Επαναφορτώσιμο / Μονό Πλαίσιο Uni-frame® Többszór használatos / Uni-frame® Ricaricabile / Uni-Frame® 再装填可能 / Uni-frame® Atkārtoti ielādējams / Uni-frame® Uždedams / Uni-frame® Gjenbrukbar ramme / Uni-frame® Do uzupelniania / Uni-frame® Recarregável/Uni-frame® Recarregável / Uni-frame® Reincàrcabil / Uni-frame® Dobijateľný / Uni-frame® Polnilni / Uni-frame® Armazón recargable / Uni-frame® Påfyllbar / Uni-frame® Yeniden Doldurulabilir / Uni-frame®</p>	 <p>Latched 已锁定 Zaključano Zachycen západkou Lást Gesloten Lukustatud Lukittu Enclenché Arretiert Ασφαλισμένο Rögzített Chiuso 掛かっている Fiksēts Užfiskuota Lást Zabezpieczony zatraskiem Preso Preso Inciuat So západkou Zaskočen Bloqueado Spårrad Mandali kapalı</p>	 <p>Unlatched 已解锁 Otključano Ödroyeno Åben Geopend Lukustamata Auki Déverrouillé Gelöst Απασφαλισμένο Rögzítetlen Aperto 外れている Nefiksēts Neuzfiksiskuota Ulást Niezatrasnięty Nåo preso Solto Descuiat Bez západky Odpahņjen Desbloqueado Ospårrad Mandali açık</p>
---	--	---	--

<p>Breast Support</p> 	<p>胸支架 Osłona za prsa Opěrka prsou Bryststøtte Borstondersteuning Rinnatugi Rintatuki Soutien poitrine Bruststütze Στήριγμα Στήθους Emlődötmasz Supporto mammella 胸部サポート</p>	<p>Krūšu atbalsts Krūtinės atrama Bryststøtte Podparcie piersi Supporte de Mama Apoio do peito Support piept Podpora hrude Podpora za prsi Supporte de pecho Bröststöd Göğüs Desteği</p>
---	--	--

	<p>I. Groin Piece 腹股沟件 Genitalna zaštita Triselný díl Inguenstykke Liesstuk Kubemeosa Nivuskappale Coquille Leistenstück Τηχηα Βουβιόνα Ágyékélem Pezzo glutei 鼠径部ピース Cirkšna dala Dalis kirksnies srityje</p>	<p>Lyskestykke Wspornik pachwinowy Peça para Virilha Peça para virilha Element pentru regiunea inghinală Diel na slabiny Opora za prepono Pieza de ingie Ljumskdæl Kasik Parçası</p>
---	--	--

INTENDED USE

- The intended use of this device is to immobilize the patient during radiation therapy.
- Devices are used to aid in the support and positioning of patients during an MRI.
- The device is indicated to aid in supporting, positioning, and/or immobilization of adult and pediatric patients undergoing radiation therapy of the head, brain, neck, and spine including radiosurgery and electron, photon, and proton treatments. The device is also used during image acquisition to support treatment planning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

MRI SAFETY INFORMATION

MR Safe

- Devices are MR Safe.
- HipFix® Thermoplastic only: In non-clinical testing, the image artifact caused by the device extends up to 12.1mm from groin piece of the HipFix® thermoplastic when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

⚠ WARNING

- Do not use if device appears damaged.
- Device should be used for a single patient throughout setup and treatment cycle.
- Ensure thermoplastic is properly aligned over patient before forming.
- Ensure nose mouth portal is aligned during molding to allow patient to breathe from nose and mouth during future use.
- Ensure patient is comfortable during molding process to ensure fit is comfortable during future use.
- Form thermoplastic to avoid contact with tracheotomy tube.
- Cutting thermoplastic masks may reduce the rigidity of the device. If cutting is necessary, ensure rigidity of mask is maintained by adding thermoplastic strips. Thermoplastic strips are not included with masks but can be purchased by contacting your CIVCO Radiotherapy representative.
- Verify all angles of treatment, attenuation characteristics and WET values prior to treating patients. Refer to Technical Data Sheet at WWW.CIVCORT.COM for WET values.
- Ensure thermoplastic is properly attached using clamps, pins or quick fasteners.
- Accessory carbon fiber frames and profiles are not MR safe.
- Thermoplastic material may cause skin irritation in individuals with sensitive skin.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

Thermoplastic Type	Waterbath Temperature	Waterbath Heating Time	Dry Heat Oven Temperature*	Dry Heat Oven Heating Time*
ZENTEC™	75°C (165°F)	2–3 min	Not approved for oven	Not approved for oven
Head Only; 25cm, 35cm Thermoplastic Strips	75°C (165°F)	4 min	75°C (165°F)	4–7 min
IMRT Reinforced Thermoplastics™ Head Only	75°C (165°F)	5 min	75°C (165°F)	8–11 min
Head, Neck and Shoulders; Breast Support; Sheet	75°C (165°F)	6 min	75°C (165°F)	4–8 min
IMRT Reinforced Thermoplastics Head, Neck and Shoulders; HipFix	75°C (165°F)	7–9 min	75°C (165°F)	8–11 min
Posicast® Thermoplastic	70°C (158°F)	1 min	75°C (165°F)	5-9 min
Mastercast™ Pro	70°C (158°F)	1 min	75°C (165°F)	5-9 min

* NOTE: Temperatures and times for dry heat oven may vary depending on model.

USING THERMOPLASTICS

⚠ WARNING

- Type-S™: Ensure pins are unlatched and parallel with edge of frame (refer to Type-S™/ Uni-frame® illustration).

1. Setup waterbath with 2–3" of water and heat or turn on dry heat oven (refer to table for appropriate temperature).
2. Insert thermoplastic into reloadable frame or profile (if applicable).
3. Place thermoplastic into waterbath/dry heat oven for directed time.

NOTE: Improper heating time may cause difficult application.

- Remove thermoplastic from waterbath/dry heat oven.


⚠ WARNING

- Water is hot, do not reach into water to remove device.
- Feel thermoplastic mask to make sure it is not too hot to come in contact with patient's skin.

- Quickly remove excess water from thermoplastic (*less than five seconds*).
- Position thermoplastic over treatment area.

NOTE: • Overstretching of thermoplastic may affect rigidity of mask. After pulling thermoplastic mask, ensure rigidity of mask is maintained.
 • Some thermoplastic masks may require assistance in pulling of mask.

- Gently pull sides of frame outward and down to attachment base.

- Use all available clamps, pins or quick fasteners (*as appropriate according to baseplate*) to latch thermoplastic to baseplate.

NOTE: Ensure patient's hair is not caught between mask and baseplate or extension.

- Gently mold material to contours of treatment area.
- Allow thermoplastic to cool completely before releasing from attachment base.


USING CLEARVISION™, CLEARVISION™ 2 AND PUREVISION™ THERMOPLASTIC MASKS

- Position mask over patient, aligning mask to capture forehead and chin with the thermoplastic.

NOTE: • Two people are required to pull mask, one to pull mask, while the other assists in securing the mask at forehead and chin locations.
 • While pulling mask ensure thermoplastic is not folded incorrectly under itself.


- Gently pull down to appropriate attachment base (*refer to attachment base instruction*); ensuring thermoplastic is secure at forehead and chin locations.

NOTE: Ensure desired anatomical locations are free of thermoplastic material as required by camera system region of interest.


USING PUREVISION ONLY (STEPS 3-5):

3. Fold side flaps and chin area down and away from face to create an open mask.

4. Center 35cm thermoplastic strip at medial point of chin and wrap strip superiorly toward base end of skull, joining ends together at back of mask. Keep strip below cheekbones for optimal openness.
5. Center 25cm thermoplastic strip across forehead and pull vertically downward along lateral sides of mask, making sure the strip is embedded into mask.


- I. Reloadable Mask
- II. 35cm Thermoplastic Strip
- III. 25cm Thermoplastic Strip
- IV. Frame

RELEASING QUICKFASTENERS FOR POSICAST®

1. Press profile in the middle with one hand. Maintain this pressure through steps 1-4.
2. Push first quickfastener with other hand.
3. Pull first quickfastener "CLICK" should be audible.
4. Push and pull second quickfastener.
5. Release pressure from profile. Pull the profile from the baseplate in vertical direction.
6. Release all profiles one by one as described above.

NOTE: If it is difficult to remove the quickfasteners, use a toothpick or a cotton swab and carefully remove any residue in the quickfastener holes.

NOTE: To reorder components listed, contact your CIVCO Radiotherapy sales representative.

- Thermoplastic, Reloadable (760030, 760145, 770145, 910030, 910044, 910145, 920033, 920145, 960030, 960145)

REPROCESSING

⚠ WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. If necessary, clean surface of thermoplastic and/or frames and pins by removing gross visual contaminants with alcohol or warm soapy water. If gross visual contamination cannot be removed repeat cleaning steps and if necessary, discard device.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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800.842.8688, +1 712.737.8688

info@civcort.com

COPYRIGHT © 2023 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. CIVCO RADIOTHERAPY AND CIVCO RT ARE TRADEMARKS OF MEDTEC LLC. HIPFIX AND UNI-FRAME ARE REGISTERED TRADEMARKS OF MEDTEC LLC. CLEARVISION, IMRT REINFORCED THERMOPLASTIC, MASTERCAS, PUREVISION, TYPE-S AND ZENTEC ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.

WWW.CIVCORT.COM