
C-Qual™


REFERENCE GUIDE

CIVCO
Radiotherapy™

CE

- I. Lower Board
- II. Bottom Stop
- III. Tilting Board
- IV. Head Support
- V. Upper Arm Support
- VI. Lower Arm Support
- VII. Lower Arm Support Knob
- VIII. Sure-Set Angulation - Low Angle Support
- IX. Lock
- X. Sure-Set Angulation - High Angle Support
- XI. Thermoplastic Mask Attachment Points
- XII. Upper Arm Support Lower Knob
- XIII. Upper Arm Support Upper Knob
- XIV. Ruler

INTENDED USE

Device is intended to achieve a desired set-up for a breast-patient for radiotherapy and diagnostics.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

- Do not use if device appears damaged.
- Do not use device with MRI scanner.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- Ensure all accessories are attached securely prior to treatment.
- Do not allow patient to reposition themselves.
- Ensure patient remains still for duration of setup and treatment.
- Do not use Lok-Bar™ to secure positioning devices.
- When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at WWW.CIVCORT.COM.
- Verify patient position with completed setup sheet prior to treatment.
- Do not exceed patient weight of 330 lbs. (150 kg).
- To prevent damage, grasp both the tilting board and lower board when handling.
- Do not reposition device with patient on it.
- Verification image should be taken to confirm patient position prior to treatment.
- Only use CIVCO Radiotherapy compatible accessories.
- This product is manufactured with carbon fiber, and rough handling may cause damage resulting in fibers coming off. If this happens, stop using the product to prevent fibers from getting in contact with the patient's or user's skin.


NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

ATTACHING C-QUAL™

⚠ WARNING

- Lok-Bar™ is a positioning device only, do not use to secure patient.


1. Snap Lok-Bar™ at appropriate longitudinal position. Ensure Lok-Bar™ is secure.
2. Place C-Qual™ onto Lok-Bar™ pins carefully to prevent damage. Ensure C-Qual™ and Lok-Bar™ are secure.

ADJUSTING C-QUAL™

NOTE: When positioning patient for the first time, use C-Qual™ setup sheet to record all adjustments.



1. Pull lock to unlock tilting board from lower board.


2. Position tilting board:

WARNING


- Ensure tilting board is fixed properly before placing patient.

 <p>I.</p> <p>For 5°, 7.5°, 10°, 12.5°, 15° angles lift low angle support and secure tilting board into appropriate location.</p>	 <p>II.</p> <p>For 17.5°, 20°, 22.5°, 25° angles lay low angle support onto lower board. Separate high angle support from tilting board and select appropriate location along lower board.</p>	<p>I. Sure-Set Angulation - Low Angle Support II. Sure-Set Angulation - High Angle Support</p>
--	--	---




3. Position head support as needed (H1-H5).


4. Slide or lift on bottom stop to achieve appropriate location (B1-B13).

5. Adjust upper arm support position:

 <p>Loosen upper arm support lower knob to rotate arm support to desired location (A-P).</p>	 <p>Slide upper arm support for appropriate horizontal position (T-Z). Tighten upper arm support lower knob.</p>	 <p>Loosen upper arm support upper knob. Adjust angle of arm support by moving knob along curved scale to appropriate position (1-12). Tighten knob.</p>
---	---	---

6. Loosen lower arm support knob to adjust lower arm support height (A1-A2). Tighten knob.

7. Head and breast thermoplastics can be attached to C-Qual™ by inserting into attachment points (See thermoplastics instruction).
 8. Prior to removing C-Qual™ from couchtop, lock tilting board to lower board.

REPROCESSING

⚠ WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.
- Inspect device for damage such as cracks or breakage. If damage is evident, discontinue use and contact CIVCO Radiotherapy for further instructions as needed.


1. Wipe all surfaces with common germicidal or antiseptic wipe, such as alcohol. If necessary, disinfect with a mild, non-abrasive cleanser.
 2. Wipe dry with lint-free cloth or allow to air dry prior to reassembly and use.


MAINTENANCE

- NOTE:
- Inspect device periodically for signs of damage and general wear.
 - Inspect arm supports for damage and wear. Clean to remove grease and contaminants from arm support components.

STORAGE

- NOTE:
- Wrap device in clean soft cloth to prevent scratching. Store in cool, dry place.
 - Before C-Qual™ storage holder installation, consult your facility's manager.
 - Device should be stored on C-Qual™ storage rack.

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CIVCO

Radiotherapy™

 CIVCO RADIOTHERAPY
1401 8th Street SE
Orange City, IA 51041
United States
info@civcort.com

COPYRIGHT © 2023 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. CIVCO RADIOTHERAPY AND CIVCO RT ARE TRADEMARKS OF MEDTEC LLC. C-QUAL, LOK-BAR AND PRODIGY ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.

WWW.CIVCORT.COM