
# Now Validated on TrueBeam™ and VitalBeam™

CIVCO Radiotherapy is pleased to announce that Universal Couchtop™ Long Extension (LE) is now validated and verified on Varian's TrueBeam and VitalBeam™ Radiotherapy Systems.

Validated on TrueBeam version 2.7 and higher software, including PerfectPitch™ 6 DoF and 4 DoF couch.

Please contact your Varian sales representative for specific configuration.



# Reduce set-up and treatment errors with dosimetrically matched Universal Couchtops



Available for most major CT and Linac models

## Seamless transition from simulation to treatment

Installing Universal Couchtop on all department pedestals will offer the same:

- Dose attenuation
- Indexing
- Deflection
- · Immobilization device attachment



# **Extensions to Treat Many Disease Sites**



Body Pro-Lok ONE™ O Body Pro-Lok™

All Long Extensions (SBRT Immobilization with Bridge or Belt option)



trUpoint I ARCH™ Posifix Extension

(SRS/SRT Immobilization)



Type-S™ & Posifix® Extensions (Head & Neck Positioning)



Universal Couchtop Long Extensions have been validated on TrueBeam™ and VitalBeam™. CIVCO positioning and immobilization devices have not been validated with TrueBeam or VitalBeam and Varian is not liable for use of these products. Universal Couchtop has not been validated on HyperArc™ or Edge™.

## **UCT LE Statistics**

#### **Attenuation** 6MV

Gantry	Transmission (10 x 10)	Transmission (5 x 5)
180°	0.971	0.969
200° / 160°	0.968	n/a
210°	0.967	n/a
220° / 140°	0.962	n/a

### **Load Rating**

Entire System: 220 kg/484 lbs Extension Only: 178 kg/392 lbs

Attenuation and surface dose values are approximate due to variation in machine characteristics.

#### **Surface Dose**

	6MV 10 x 10	6MV 20 x 20	18MV 10 x 10	18 MV 20 x 20
Open Field	19.1%	30.5%	16.1%	20.6%
UCT LE	94.1%	97.2%	69.9%	79.4%
Posi Extension P1	82.6%	89.7%	56.9%	69.4%
Posi Extension P2	89.6%	93.2%	64.6%	73.6%



Global Sales Office

Orange City, Iowa, USA 800.842.8688 | +1 712.737.8688 info@CivcoRT.com www.CivcoRT.com

