

Body Pro-Lok ONE™ ONEPlatform™

TECHNICAL DATA SHEET




PRODUCT NAME	ONEPlatform
PRODUCT NUMBER	BPL003
GENERAL DESCRIPTION	An integral part of the Body Pro-Lok ONE System, the ONEPlatform features the same profile and indexing options as the Universal Couchtop™, providing seamless integration of positioning devices via three-pin Lok-Bar™ indexing.
INDICATIONS FOR USE	The device is intended to aid in support and positioning during MR, radiological and other procedures.
MR SAFE	Yes
INTEGRATION	2 or 3-Pin indexable to couchtop
FDA CLEARANCE	Yes
CE MARKED (EU DOC)	Yes
OTHER COUNTRY REGISTRATIONS	Available upon request
COUNTRY OF MANUFACTURE	USA

TECHNICAL DETAILS

COMPONENT	BPL003 ONEPlatform	
PRIMARY MATERIALS	Kevlar, Foam Core, Phenolic Inserts	
WEIGHT (APPROX)	7.5 kg / 16.5 lbs	
MAX PATIENT WEIGHT	226.8 kg / 500 lbs	
DIMENSIONS (CM)	182.9 L x 53 W x 3.5 H (max)	
DOSIMETRIC PROPERTIES	LOCATION	ATTENUATION 6 MV PHOTONS 10 X 10 FIELD SIZE (CM ²)
	Center of platform	2.1%
	Performance Specifications are for reference only, values may vary based on setup and equipment.	

TECHNICAL DRAWING



PACKAGING

UNIT PACK	1 Piece
UNITS / BOX	1 Piece
SHELF LIFE	Not Applicable
WARRANTY	See website for Product Warranty Statement
STERILIZATION	Non-Sterile
RELEASE DATE	September 2014

ACCESSORIES

SYSTEM CUSTOMIZATION	Body Pro-Lok ONEBridge (BPL001) Body Pro-Lok ONE Respiratory Belt (BPL026) Additional devices compatible with 2 and 3-pin Lok-Bars/Locating Bars
----------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------