

# Chabner XRT® Radiation Bra

## CIVCO RT Technical Data Sheet



<b>Product Name</b>	Chabner XRT® Radiation Bra			
<b>Product Number</b>	<b>Individual Size</b>	<b>Radiation Bra Qty 1</b>	<b>Radiation Bra Qty 5</b>	
	1	CRTBRA01	CRTBRA01PK5	
	2	CRTBRA02	CRTBRA02PK5	
	3	CRTBRA03	CRTBRA03PK5	
	4	CRTBRA04	CRTBRA04PK5	
	5	CRTBRA05	CRTBRA05PK5	
	6	CRTBRA06	CRTBRA06PK5	
	7	CRTBRA07	CRTBRA07PK5	
	8	CRTBRA08	CRTBRA08PK5	
	9	CRTBRA09	CRTBRA09PK5	
	10	CRTBRA10	CRTBRA10PK5	
	11	CRTBRA11	CRTBRA11PK5	
	12	CRTBRA12	CRTBRA12PK5	
	13	CRTBRA13	CRTBRA13PK5	
<b>Kits that include this product</b>	CRTBRASPK1-5 CRTBRAKIT1-5	CRTBRASPK5-9 CRTBRAKIT5-9	CRTBRASPK9-13 CRTBRAKIT9-13	CRTBRASPK1-13 CRTBRAKIT1-13
<b>General Description</b>	Provides optimal breast support during simulation, planning and radiation treatment			
<b>Indications for use</b>	This device is indicated to aid in supporting adolescent and adult patients undergoing radiation therapy including electron, photon, and proton treatments. The device is also used during image acquisition to support treatment planning. The device is not intended for use with patients under 12 years of age.			
<b>MRI Compatible</b>	Yes			
<b>FDA Clearance</b>	Yes			
<b>CE Marked (EU DOC)</b>	Yes			
<b>Other Country Registrations</b>	Available upon request			
<b>Country of Manufacture</b>	USA			

# Technical Details

## Primary Materials


UBL Tricot, Mesh Warp Knit, Synthetic Elastomer, Nylon, Polyester

## Dosimetric Properties

Material	6 MV photons	6 MeV electrons	Proton water equivalency R80
Tricot	0%	0%	–
TPU ●	0%	0%	0.2 mm
Elastic	0%	.1%	–
Velcro ●	0%	.7%	–
Mesh ●	0%	0%	0.1 mm

Performance Specifications are for reference only, values may vary based on setup and equipment.

## Technical Drawing



Part Number	Dimensions: Total Bottom Band on the Half (Fully Closed)	Product Weight (Approx)
CRTBRA01	14.25 in (36.2 cm)	2.1 oz (60 grams)
CRTBRA02	14.25 in (36.2 cm)	2.3 oz (65 grams)
CRTBRA03	15.75 in (40 cm)	2.5 oz (70 grams)
CRTBRA04	15.75 in (40 cm)	2.6 oz (75 grams)
CRTBRA05	16.00 in (40.6 cm)	2.8 oz (80 grams)
CRTBRA06	17.5 in (44.5 cm)	3.0 oz (85 grams)
CRTBRA07	17.375 in (44.1 cm)	3.2 oz (90 grams)
CRTBRA08	18.875 in (47.9 cm)	3.4 oz (95 grams)
CRTBRA09	18.75 in (47.6 cm)	3.5 oz (100 grams)
CRTBRA10	18.75 in (47.6 cm)	3.5 oz (100 grams)
CRTBRA11	18.75 in (47.6 cm)	3.7 oz (105 grams)
CRTBRA12	20.25 in (51.4 cm)	4.1 oz (115 grams)
CRTBRA13	20.25 in (51.4 cm)	4.2 oz (120 grams)

## Packaging

Unit Pack	Each
Units / Box	One device
Shelf Life	n/a
Warranty	<a href="#">See CIVCO RT website for Product Warranty Statement</a>
Sterilization	Non-Sterile
Launch Date	October 2017 – USA March 2018 – Global

### Global Sales Office

Orange City, Iowa, USA  
800.842.8688 | +1 712.737.8688

info@CivcoRT.com  
www.CivcoRT.com

