

Respiratory Plate

CIVCO RT Technical Data Sheet




Product Name	Wide Respiratory Plate	Standard Respiratory Plate
Product Number	BPL008	MTSBRT003
Kits that include this product	BPL100, BPL110, BPL120, BPL200, BPL210, BPL220, BPL230	BPLONEUPG, BPLONEUPG2, BPLONEUPG3, MTSBRT001, MTSBRT001L, MTSBRT320, MTSBRT420K, MTSBRT350
General Description	When attached to a Body Pro-Lok™ Bridge or ONEBridge™, the Standard Respiratory Plate and Wide Respiratory Plate place pressure at the level of the diaphragm to assist in restricting respiratory movement the influences tumor movement.	
Indications for use	Lung and Abdominal Stereotactic Radiation Body Therapy (SBRT)	
MRI Compatible	Yes	
Integration	Body Pro-Lok Bridges and ONEBridge	
FDA Clearance	Yes	
CE Marked (EU DOC)	Yes	
Other Country Registrations	Available upon request	


Technical Details

	Wide Respiratory Plate	Standard Respiratory Plate
Primary Materials	Black Acetal, Black G10	Black Acetal, Black G10
Weight (approx)	0.7 lbs (0.32kg)	0.8 lbs (0.36kg)
Dimensions (mm)	258 L x 144 W x 136 H (max)	171 L x 132 W x 130 H
Country of Manufacture	USA	USA

Technical Drawing - BPL008



Technical Drawing - MTSBRT003



Packaging

	Wide Respiratory Plate	Standard Respiratory Plate
Unit Pack	1 Piece	1 Piece
Units / Box	1 Piece	1 Piece
Shelf Life	n/a	n/a
Warranty	See CIVCO RT website for Product Warranty Statement	
Sterilization	Non-sterile	Non-sterile
Launch Date	December 2016	2008

Accessories

	Wide Respiratory Plate	Standard Respiratory Plate
System Customization	ComfortCare™ Cushion BPL007	n/a

Global Sales Office

Orange City, Iowa, USA
800.842.8688 | +1 712.737.8688

info@CivcoRT.com
www.CivcoRT.com

