GrayDuck Dart™

ORAL POSITIONING STENT

NEW & IMPROVED DESIGN



TONGUE-DISPLACING AND IMMOBILIZING STENT

GrayDuck Dart is currently the only custom fitted bite block system that can position the tongue in one of several positions; laterally (left or right), depressing (downward) or a combination thereof.

- Tongue-displacing and immobilizing oral stent
- · Capable of providing repeatable tongue deviation and immobilizing in several directions
- · Designed to provide repeatable immobilization and tissue sparing during radiotherapy
- Ergonomically developed to be a quick, easy-to-fit custom oral stent

SMALLER AND EASIER-TO-USE DESIGN

The GrayDuck Dart offers a significantly smaller and easier-to-use stent while retaining the same level of customization as the original GrayDuck Stent™ design. Additional benefits include:

- Smaller stent size increases patient tolerance and improved insertion for patients with limited bite opening
- Three simple steps allow for faster workflow
- · Spacers allow for increased bite opening, if needed
- · Universal stent orientation allows placement in patient's mouth on either arch tray side



GrayDuck Dart is 55% smaller than the original stent design



GrayDuck Dart. Each Kit includes dart stent, 2 paddles, and 2 spacers.

GDS003

GrayDuck Dart, Qty 4

CE MARKED, 510(K) CLEARED. US PATENT #9,504,537 AND OTHER PATENTS PENDING.



