

# Universal Couchtop™ for Safe, Efficient Treatment Planning & Delivery

 ProtonSeries™

## Junction, Rail and Artifact-free

Specifically designed to answer the demands of Proton, the Universal Couchtop Long Extension contains no rails or junctions in the treatment area to help facilitate better dose delivery and target coverage.

### Universal Couchtop™ Long Extension

- Large treatment area of 58.9" (149.7 cm)
- Homogenous design for consistent treatment planning & delivery
- Utilizes Prodigy™ 2 indexing at 7 cm increments
- Low attenuation and surface dose
- Self-locking extension seam located outside of treatment field
- 5 Year Warranty



### IBA Systems Proteus® ONE and Proteus PLUS

MTIL6840 F-Series

MTIL6841 L-Series

Both part numbers include: UCT LE Couchtop with Rectangular Extension (MTIL6673) and Prodigy™ 2 Lok-Bar™



Global Sales Office  
Orange City, Iowa, USA  
800.842.8688 | +1 712.737.8688

info@CivcoRT.com  
www.CivcoRT.com



**CIVCO**®  
Radiotherapy