## Type-S™ Lok-Bar™

## for use with Type-S Overlays or Extensions, including TotalRT All in One Platform

### CIVCO RT Technical Data Sheet




Product Name	Type-S Lok-Bar
Product Number	MTWBL001
General Description	The Type-S Lok-Bar allows Wing Boards™, Monarch™ Overhead Arm Positioner and other positioning and immobilization devices to be indexed in the head & neck area on a Type-S Overlay or Extension, including the TotalRT All in One Platform.
Indications for Use	The device is intended to aid in supporting and positioning adult and pediatric patients undergoing radiation therapy including electron, photon, and proton treatments. The device is also used during image acquisition to support treatment planning.
MRI Compatible	MR Safe
Integration	Place in shoulder holes of Type-S overlay or extension, including TotalRT All in One (TRT100)
FDA Clearance	Yes
CE Marked (EU DOC)	Yes
Other Country Registrations	Available upon request
Country of Manufacture	USA




## Technical Details

Primary Materials	Black Acetal
Weight (approx)	.06 kg / .13 lbs
Dimensions (mm)	431.8 L x 25.40 W x 25.46 H (max)

#### **Technical Drawing - MTWBL001**







# Packaging

Unit Pack	Each
Units/Box	One device, may have several component parts in box
Shelf Life	n/a
Warranty	See CIVCO RT website for Product Warranty Statement
Sterilization	Non-Sterile
Release Date	October 2021



Orange City, Iowa, USA 800.842.8688 | +1 712.737.8688 info@CivcoRT.com www.CivcoRT.com

