

# AccuForm™ Shim

## CIVCO Technical Data Sheet



<b>Product Name</b>	AccuForm™ Shim
<b>Product Number</b>	MTACLSHIM
<b>General Description</b>	The Shim is an accessory for use after forming the patient's AccuForm cushion. Using this shim can counteract patient swelling and/or suspected mask shrinkage, increasing patient comfort while maintaining ideal immobilization.
<b>Indications for use</b>	The device is indicated to aid in supporting, positioning, and/or immobilization of adult and pediatric patients
<b>Contraindications</b>	None
<b>MRI Compatible</b>	MR Safe
<b>FDA Clearance</b>	Yes
<b>CE Marked (EU DOC)</b>	Yes
<b>Other Country Registrations</b>	Available upon request

---

# Technical Details

---

**Primary Materials** Corss-linked Polyolefin Foam

---


**Weight (approx)** 0.009 kg/.02 lbs

---

**Dimensions (mm)** 229 L x 93 W x 1.57 H (min.)

---

## Technical Drawing



---

# Packaging

---

<b>Unit Pack</b>	10
<b>Units / Box</b>	Qty 1
<b>Shelf Life</b>	n/a
<b>Warranty</b>	<a href="#">See CIVCO website for Product Warranty Statement</a>
<b>Sterilization</b>	Non-Sterile
<b>Launch Date</b>	March 2021

---

## Global Sales Office

Orange City, Iowa, USA  
800.842.8688 | +1 712.737.8688

[info@CivcoRT.com](mailto:info@CivcoRT.com)  
[www.CivcoRT.com](http://www.CivcoRT.com)

