

Increased Patient Comfort and Rigidity

Solstat™
Immobilization System

Solstat™ Immobilization System

CIVCO Radiotherapy brings leading-edge innovation, increased patient comfort and rigidity for head and neck treatment with the Solstat Immobilization System.

The Solstat system is comprised of a rigid head support, customizable cushion or non-custom headrest and a dedicated thermoplastic mask designed for secure and simple attachment. These components create a highly customizable and comfortable setup for more accurate head and neck treatments.


Not available in the U.S.



Featuring IMRT
Reinforced™ Thermoplastic



Features & Clinical Benefits

- Midline thermoplastic mask attachment provides increased rigidity and minimizes shrinkage to enhance patient comfort
- Highly conformal to the patient when using a customizable cushion and dedicated thermoplastic mask
- Adaptable with Precise Bite to aid in effective and repeatable positioning
- Unique thermoplastic frame provides easy to use, secure and quiet mask attachment
- Type-S™ and Posifix® version compatible
- MR Compatible  **MRSeries™**



Featuring ClearVision™
Open-face Thermoplastic with
Precise Bite™, compatible with
VisionRT's AlignRT® camera system



See CIVCO's YouTube channel for
the Solstat Instructional Video

System Customization

The system is comprised of a rigid head support, customizable cushion or non-custom headrest and a dedicated thermoplastic mask designed for secure and simple attachment. These components create a highly customizable, comfortable and rigid setup for accurate head and neck treatments.

Individual Order Items

Head Support

CHS05	Posifix® (left image)
CHS06	Type-S™ (right image)



ClearVision™ Thermoplastic

visiφnrt | COMPATIBLE

MTAPCID1832SG	Standard (left)
MTAPCIN1832SG	Narrow (right)




IMRT Reinforced Thermoplastic™, Style 18

MTAPCID1832	Standard (left)
MTAPCIN1832	Narrow (right)



Headrests

CHS40K	Headrest Kit (includes B, C, & F), Qty 3
CHS40	Headrest B, Black, Qty 1
CHS41	Headrest C, Black, Qty 1
CHS42	Headrest F, Black, Qty 1




AccuForm™ Cushions

MTACL2035HF	20 x 35 cm Chamfered, Half Fill
MTACL2035	20 x 35 cm Chamfered



Precise Bite™

MTCBPBITE	Standard, Qty 5
MTCBPBITES	Small, Qty 5



- Vision RT is unable to comment on, nor make claims, about the immobilization rigidity or patient alignment reproducibility when using this mask;
- Compatibility is subject to users following CIVCO guidelines for forming open masks in conjunction with Surface Guided Radiotherapy Systems, and subject to users following Vision RT's guidelines for drawing Regions of Interest for use in SRS with Open Masks;
- The mask vendor is responsible for ensuring that users are in receipt of guidelines for forming masks for use in conjunction with Surface Guided Radiotherapy Systems. Vendors may change their guidelines from time to time without reviewing these changes with Vision RT. Vision RT is not liable for misuse of masks;
- Vision RT makes no comment on, and is not liable for, suitability or otherwise of masks for SRS or any other procedure;
- Purchasers are responsible for making their own choices as to the masks selected by them.

Global Sales Office

Orange City, Iowa, USA
800.842.8688 | +1 712.737.8688

info@CivcoRT.com
www.CivcoRT.com

CIVCO[®]
Radiotherapy