


TotalRT All in One Solution





 **TotalRT™**
All in One Solution

CIVCO
Radiotherapy™
a better way

Increase Efficiency, Position Devices Easily and Gain Confidence in Patient Setup












MRSeries products are intended for use in the MR environment. Always refer to product labeling to identify whether a device is MR Safe  or MR Conditional. 

Base Package

Base Package items are sold separately. Suggested quantities are listed. Consult your sales representative for more information.

 All Base Package devices are MR Safe.

	TotalRT Platform 201.3 cm, Glass Fiber TRT100 Qty 1
	Pro Wing Board™ Includes Pro Handles and Gel Head Rest OAP220 Qty 1
	Customizable Head Support (Type-S™) CHS02 Qty 1
	AccuForm™ Cushion 20 x 35, Chamfered MTACL2035 Qty 1
	Kneefix™ 3 301052 Qty 1
	Feetfix™ 3 301071 Qty 1
	Prodigy™ 2 MR Three-Pin Lok-Bar™ MTIL6054 Qty 3
	Type-S™ Lok-Bar MTWBL001 Qty 1
	Handles PRO360R and PRO360L Qty 1 Each

Body Pro-Lok ONE™ SBRT Immobilization

Body Pro-Lok ONE features the ONEBridge™, designed with unique lateral & tilting offset capability with variable height adjustment. A Respiratory Belt option is also available.



Body Pro-Lok ONE™







SBRT Package



Additional SBRT Packages are available. Please visit CivcoRT.com for additional information.



BPL410 Includes:



- 




TotalRT Platform
201.3 cm, Glass Fiber
Qty 1
- 



Body Pro-Lok ONEBridge™
Qty 2
- 



Wide Respiratory Plate with ComfortCare™ Cushion Starter Kit
Includes Wide Plate, 3 Cushions and Hand Pump
Qty 1
- 


Clam-Lok™ Cushion
Qty 2
- 


Kneefix™ 3
Qty 1
- 


Feetfix™ 3
Qty 1
- 


Prodigy™ 2 MR Three-Pin Lok-Bar™
Qty 3
- 


Type-S™ Lok-Bar
Qty 1
- 


ONEBridge Storage Base
Qty 2

Solstice™ SRS Package

To complete this setup, refer to the **Base Package** on page 3.
Consult your sales representative for more information.

MR All Solstice SRS Package devices are MR Safe.

CHS04MRK Includes:



Glass Fiber Head Support
Qty 1



Solstice ClearVision™ Open-face Thermoplastic Standard
Qty 10



AccuForm™ Cushion
20 x 35 cm Chamfered
Qty 10



Precise Bite™
Qty 5



TotalRT All in One Solution

Features & Benefits

- Easily indexable on CT, MR and photon treatment tables with two Lok-Bars™
- Large, artifact-free and homogeneous treatment area for precise dose delivery, including hypofractionation
- Repeatable positioning of CIVCO RT devices, even between treatment rooms
- A single therapist can carry and index to couchtop
- Wall storage device allows effective way to store and protect the platform


Specifications

- Part number: TRT100
- Lightweight: 7.3 kg / 16 lbs
- Dimensions (cm): 201.3 (L) x 61.9 (W) x 3.5 (H)
- Platform is MR Safe

MR



View the TotalRT
Technical Data Sheet!





CIVCO

Radiotherapy™



Global Sales Office

1401 8th Street SE

Orange City, IA 51041, USA

800.842.8688 | +1 712.737.8688

info@CivcoRT.com | www.CivcoRT.com

COPYRIGHT © 2023. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. HIPFIX IS A REGISTERED TRADEMARK OF MEDTEC LLC. BODY PRO-LOK, BODY PRO-LOK ONE, CIVCO RADIO THERAPY, CIVCO RT, CLAM-LOK, COMFORTCARE, FEETFIX, KNEEFIX, LOK-BAR, MRSERIES, ONEBRIDGE, PRODIGY, SOLSTICE, TYPE-S, UNIVERSAL COUCHTOP AND VAC-LOK ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.

2022P1600US REV. B