

Body Pro-Lok ONE™ SBRT Immobilization

CIVCO's Body Pro-Lok ONE™ provides one simple solution for complex SBRT patient setups. This system features the ONEBridge™, with unique design features allowing for improved patient positioning and compression plate restriction.

Advancing SBRT Respiratory Restriction



Customize the patient's treatment position and comfort level with the only SBRT bridge that features lateral and tilting offset with variable height adjustment capability.



NEW! Bringing Compassion to Compression: Wide Respiratory Plate & ComfortCare™ Cushion



CIVCO's new Wide Respiratory Plate helps diffuse pressure across the abdominal region, resulting in a more tolerable compression for the patient.

ComfortCare Cushions provide an added dimension of customized comfort during compression. Each ComfortCare Cushion features innovative and patented technology with soft and moveable polystyrene pearls combined with adjustable air chambers to reduce the amount of patient movement. The cushion forms to the contours of the patient during inflation, providing an additional level of comfort and contoured pressure between the patient and the respiratory plate.

Contact your CIVCO account manager or distributor for more information

Global Sales Office
Orange City, Iowa, USA
800.842.8688 | +1 712.737.8688

info@civco.com
www.civco.com

COPYRIGHT © 2015. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. BODY PRO-LOK ONE AND ONEBRIDGE ARE TRADEMARKS OF CIVCO. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW. 2015PI041 REV B

CIVCO[®]
Radiotherapy