Solstice™ SRS Package

To complete this setup, refer to the **Base Package** on page 3. Consult your sales representative for more information.



MR All Solstice SRS Package devices are MR Safe.

CHS04MRK Includes:



Glass Fiber Head Support

Qty 1



Solstice ClearVision™ Open-face Thermoplastic Standard

Qty 10



AccuForm™ Cushion 20 x 35 cm Chamfered

Qty 10



Precise Bite™

Qty 5



TotalRT All in One Solution

Features & Benefits

- Easily indexable on CT, MR and photon treatment tables with two Lok-Bars™
- Large, artifact-free and homogeneous treatment area for precise dose delivery, including hypofractionation
- Repeatable positioning of CIVCO RT devices, even between treatment rooms
- A single therapist can carry and index to couchtop
- Wall storage device allows effective way to store and protect the platform

Specifications

- Part number: TRT100
- Lightweight: 7.3 kg / 16 lbs
- Dimensions (cm): 201.3 (L) x 61.9 (W) x 3.5 (H)
- Platform is MR Safe





View the TotalRT Technical Data Sheet!

