

# CIVCO RT Set-up Sheet

## MT-350-N, -W Breastboard


Patient: \_\_\_\_\_ Date: \_\_\_\_\_

Notes:


1. The Lok-Bar™ is indexed to the table at (i.e. H1, F2): \_\_\_\_\_

### 2. Elevation (MT-350-N or -W only)

| Pin# | Bent | Straight |
|------|------|----------|
| 1    | 5°   | X        |
| 2    | 7.5° | 12.5°    |
| 3    | X    | 15°      |
| 4    | X    | 17.5°    |
| 5    | 10°  | 20°      |
| 6    | X    | 25°      |




### 4. Arm Rotation




Right # \_\_\_\_\_ Inner Outer  
 Left # \_\_\_\_\_ Inner Outer

### 3. Head Cup



A B C D E F  
 G H I J K L  
 M N O P Q

### 5. Upper Arm Length




Right: A B C D E

Left: A B C D E

See next page


# Set-up Sheet - Breastboard, MT-350-N, -W

## 6. Upper Arm Incline




Right # \_\_\_\_\_ Left # \_\_\_\_\_

## 9. Wrist Support Locking Clamp



Right 1 2 3 4 5 6 7 8 9 10 11  
 Left 1 2 3 4 5 6 7 8 9 10 11

## 7. Wrist Support Placement




Right 1 2 3 4 5 Left 1 2 3 4 5

## 10. Hip Stop




A B C D E F G H I J K L M N O

## 8. Wrist Support Adjustable Post



Right 1 2 3 4 5 6 7  
 Left 1 2 3 4 5 6 7



Right A B C D  
 Left A B C D

Global Sales Office  
 Orange City, Iowa, USA  
 800.842.8688 | +1 712.737.8688

info@CivcoRT.com  
 www.CivcoRT.com

