

C-Qual M™ Breastboard and/or Monarch™ Overhead Arm Positioner

Patient:

Date:


Notes:

Head Support Used


1 Index Lok-Bar™ to table (ex: H1, F2)

If using Monarch with Lok-Bar only, skip to Step 5

2 Set Tilting Board Elevation Angle (5° - 25°)



3 Position Bottom Stop (B0 - B10)




Contact CIVCO RT at info@CivcoRT.com for assistance in customizing this setup sheet for your clinical use.

Global Sales Office

Orange City, Iowa, USA
800.842.8688 | +1 712.737.8688

info@CivcoRT.com
www.CivcoRT.com

4 Set Monarch Overhead Arm Positioner (P0 - P77)



5 Hand Grip Used (check one)

One Piece Placement (G1 - G12)	Hand Position Top Side
-----------------------------------	--------------------------------

Individual Grips

Left (G1 - G12):

Right (G1 - G12):

6 Laser Alignment Location

7 Thermoplastics Used

Head-Only	Yes	No
Breast	Yes	No